

low form beakers glass



- ▶ Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- ▶ They meet the highest quality standards according to ISO 3819 and DIN 12331.
- ▶ These beakers have thick walls and heavy bases.
- ▶ Reinforced rims and spouts to prevent breakage during washing and transportation.
- ▶ Inscriptions are in white enamel which is specially designed for glass.

glass beakers white scale	volume	approximate height	Diameter	quantity in box
025.01.025	25 ml	50 mm	34 mm	10 pieces
025.01.050	50 ml	60 mm	42 mm	10 pieces
025.01.100	100 ml	70 mm	50 mm	10 pieces
025.01.150	150 ml	80 mm	60 mm	10 pieces
025.01.250	250 ml	95 mm	70 mm	10 pieces
025.01.400	400 ml	110 mm	80 mm	10 pieces
025.01.600	600 ml	125 mm	90 mm	10 pieces
025.01.800	800 ml	135 mm	100 mm	10 pieces
025.01.901	1000 ml	145 mm	105 mm	10 pieces
025.01.902	2000 ml	185 mm	130 mm	6 pieces
025.01.905	5000 ml	270 mm	170 mm	1 piece
025.01.910	10000 ml	325 mm	228 mm	1 piece

low form beakers polypropylene



- ▶ Manufactured from ultra clear virgin polypropylene which provides thick walls with superior chemical resistance.
- ▶ They meet the highest quality standards according to ISO 7056.
- ▶ Inscriptions are either blue enamel which is specially designed for plastic or with embossed scale.

P.P beakers embossed	P.P beakers blue scale	volume	diameter	height	quantity in pack
026.01.025	026.03.025	25 ml	32 mm	47 mm	10 pieces
026.01.050	026.03.050	50 ml	39 mm	60 mm	10 pieces
026.01.100	026.03.100	100 ml	50 mm	70 mm	10 pieces
026.01.250	026.03.250	250 ml	67 mm	93 mm	10 pieces
026.01.500	026.03.500	500 ml	88 mm	119 mm	6 pieces
026.01.901	026.03.901	1000 ml	109 mm	147 mm	6 pieces
026.01.902	026.03.902	2000 ml	133 mm	184 mm	4 pieces